



## Power Hi-Lo, Power Back Exam Table with Foot Control Model 4040-650-100

### TABLE SPECIFICATIONS

**Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat  
Max. 37" at seat (93.9cm)  
**Width:** 27" (68.6cm)  
**Length:** 57" – 76" with foot section extended (144.7cm - 193cm)  
**Leg Rest:** 21" W x 13" D (53.3cm W x 33cm D)  
**Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)  
**Stirrups:** Four Lateral Positions, Adjustable Length  
**Dual Electro-Mechanical Actuator System** for Efficient Speed  
**Electrical:** 120 VAC, 60HZ  
**International Model 4040-650-101:** 220/240V, 50/60Hz

### ACCESSORIES

042 –Multi-Positioning Patient Assist Armrests/Safety Rails  
Pair (Field Installable)  
251 –Bierhoff knee crutches  
467 –Mounting Plate for Welch Allen Light # 48955 (Patient Right)  
468 –Mounting Plate for Welch Allen Light # 48955 (Patient Left)  
822 –Stainless Steel Drain Pan with Spout (Factory installed only)  
870 –Drawer Warmer 110V (Factory Installed only)  
872 –Drawer Warmer 220V (International Models, Factory Installed only)  
900 –Lifting bracket kit  
920 –Easy installation kit  
Fire retardant foam (CAL133) available

### KEY FEATURES

- Power Hi-Lo + Power Backrest
- 650lb (294.84kg) Patient Weight Capacity
- Low Access ADA Height 18.5" ( ±0.5") 46.9cm (±1.27cm)
- Two-function Foot Control – Power Hi-Lo, Power Back
- Removable Seamless Upholstered Cushion
- Large Front Storage Drawer with Comfort Leg Pad
- Sturdy Steel Extension Leg Pullout with Seamless, Ultra-Comfort Leg Cushion
- Stirrups
- Pelvic Tilt (7°)
- Drain Pan
- Hospital Grade Electric Receptacle with Child Safety Cover (Patient Left)
- 3 Year Limited Warranty

### COLORS

