



Hand Control

**Power Hi-Lo, Power Back Exam Table with WheelBase[®], Programmable Hand Control, Foot Control
Model 4040-650-300**

TABLE SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
Width: 28" (71.12cm)
Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
Width: 28" (71.12cm)
Length: 57" – 76" with foot section extended (144.7cm - 193cm)
Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
Stirrups: Four Lateral Positions, Adjustable Length
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60 HZ
International Model 4040-650-301: 220/240V, 50/60Hz
 4040-650-301: 220/240V, 50/60Hz

KEY FEATURES

- Power Hi-Lo + Power Backrest
- Two-function Foot Control – Power Hi-Lo, Power Back
- 650lb (294.84kg) Patient Weight Capacity
- OneTouch WheelBase[®] Integrated Central-Locking Caster System
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- Hand Control with OneTouch Patient Positioning[®] 2 Programmable Positions & Return-to-Home
- Removable Seamless Upholstered Cushion
- Side Storage Drawer
- Large Front Storage Drawer with Comfort Leg Pad
- Drawer Warmer
- Sturdy Steel Extension Leg Pullout with Seamless Ultra-Comfort Leg Cushion
- Stirrups
- Pelvic Tilt (7°)
- Stainless Steel Drain Pan
- Hospital Grade Electric Receptacle with Child Safety Cover (Patient Left)
- 3 Year Limited Warranty

ACCESSORIES

- 042 Multi-Positioning Patient Assist Armrests / Safety Rails (Field Installable)
 - 251 Bierhoff Knee Crutches (Pair/Field installable)
 - 467 Mounting Plate for Welch Allen Light # 48955 (Patient Right)
 - 468 Mounting Plate for Welch Allen Light # 48955 (Patient Left)
 - 820 Stainless Steel Square Drain Pan (Factory installed only)
 - 822 Stainless Steel Drain Pan with Spout (Factory installed only)
 - 870 Drawer Warmer 110V (Factory Installed only)
 - 871 Hospital Grade Receptacle 110V (Factory installed only)
 - 900 Lifting bracket kit
- Fire retardant foam (CAL133) available

COLORS

