



Head Centric Procedure Chair with Programmable Hand Control: Model 4010-650-100

**CHAIR SPECIFICATIONS**

- Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat Max. 37" at seat (93.9cm)
- Seat Width:** 28" (71.12cm)
- Back Width:** 24.5" (62.23cm)
- Length:** 70" – 80" with foot section extended (177.8cm – 203.2cm)
- Leg Rest:** 28" W x 13" D (71.12cm W x 33.02cm D)
- Dual Electro-Mechanical Actuator System** for Efficient Speed
- Electrical:** 120 VAC, 60HZ
- International Model** 4010-650-101: 220/240V, 50/60HZ

**KEY FEATURES**

- Chair Designed for Head and Neck Centric Practices
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" ( ±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning® 2 Programmable Positions & Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

**ACCESSORIES**

- 044 Adjustable/self leveling arms pair- field installable
- 050 Removable U-shaped headrest
- 051 Removable rectangular headrest
- 467 Mounting plate for Welch Allyn Light #48955 - patient right
- 468 Mounting plate for Welch Allyn Light #48955 - patient left
- 830 Patient safety straps
- 900 Lifting bracket kit
- 920 Easy installation kit
- Fire retardant foam (CAL133) available

**COLORS**



Hand Control