



Power Hi-Lo, Power Back Exam Table with Programmable Hand Control, Foot Control: Model 4040-650-200

## TABLE SPECIFICATIONS

**Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)

**Width:** 28" (71.12cm)

**Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)

**Width:** 28" (71.12cm)

**Length:** 57" – 76" with foot section extended (144.7cm - 193cm)

**Leg Rest:** 19" W x 13" D (48.26cm W x 33cm D)

**Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)

**Stirrups:** Four Lateral Positions, Adjustable Length

**Dual Electro-Mechanical Actuator System for Efficient Speed**

**Electrical:** 120 VAC, 60 HZ

**International Model:** 220/240V, 50/60Hz

## ACCESSORIES

- 042 Multi-Positioning Patient Assist Armrests / Safety Rails (Field Installable)
  - 251 Bierhoff Knee Crutches (Pair/Field installable)
  - 467 Mounting Plate for Welch Allen Light # 48955 (Patient Right)
  - 468 Mounting Plate for Welch Allen Light # 48955 (Patient Left)
  - 820 Stainless Steel Square Drain Pan (Factory installed only)
  - 822 Stainless Steel Drain Pan with Spout (Factory installed only)
  - 870 Drawer Warmer 110V (Factory Installed only)
  - 871 Hospital Grade Receptacle 110V (Factory installed only)
  - 900 Lifting bracket kit
- Fire retardant foam (CAL133) available

## KEY FEATURES

### Mobility & Efficiencies

- Hand Control with OneTouch Patient Positioning® 2 Programmable Positions & Return-to-Home
- Two function foot control with Hi-Lo and power back
- Power Hi-Lo + Power Backrest
- Dual electro-mechanical actuator system for efficient speed

### Accessibility

- 650lb (294.84kg) Patient Weight Capacity
- Low-access ADA Height 18.5" ( ±0.5") - 46.9cm (±1.27cm)
- Pelvic tilt (7° )
- Drain Pan: Stainless steel 1.5" deep (3.81cm)
- Large front storage drawer
- Side storage drawer – reversible

### Patient Experience

- Removable, Seamless Upholstered Cushion
- Four lateral positions, adjustable length
- Comfort front leg cushion
- Heavy gauge steel extension leg pull-out to accommodate all patient heights

### Electrical

- Hospital grade receptacle 110V with child safety cover (patient left)

### Quality

- 3 year limited warranty
- Durable galvanized steel construction
- Easy-clean powder coat enamel surfaces
- International Models - 220/240V, 50/60Hz

## COLORS

