



Hand Control

Power Hi-Low, Manual Back Exam Table with Programmable Hand Control, Foot Control: Model 4070-650-200

CHAIR SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)

Width: 28" (71.12cm)

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)

Width: 28" (71.12cm)

Length: 57" – 76" with foot section extended (144.7cm - 193cm)

Leg Rest: 21" W x 13" D (53.3cm W x 33cm D)

Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)

Stirrups: Four Lateral Positions, Adjustable Length

Dual Electro-Mechanical Actuator System for Efficient Speed

Electrical: 120 VAC, 60 HZ

International Model: 220/240V, 50/60Hz

ACCESSORIES

042 Multi-Positioning Patient Assist Armrests / Safety Rails (Field Installable)

251 Bierhoff Knee Crutches (Pair/Field installable)

467 Mounting Plate for Welch Allen Light # 48955 (Patient Right)

468 Mounting Plate for Welch Allen Light # 48955 (Patient Left)

820 Stainless Steel Square Drain Pan (Factory installed only)

822 Stainless Steel Drain Pan with Spout (Factory installed only)

870 Drawer Warmer 110V (Factory Installed only)

871 Hospital Grade Receptacle 110V (Factory installed only)

900 Lifting bracket kit

Fire retardant foam (CAL133) available

KEY FEATURES

Mobility & Efficiencies

- Hand Control with OneTouch Patient Positioning® 2 Programmable Positions & Return-to-Home
- Power Hi-Lo with smooth moving articulating back adjustment
- Dual electro-mechanical actuator system for efficient speed

Accessibility

- 650lb (294.84kg) Patient Weight Capacity
- Low-access ADA Height 18.5" (±0.5") - 46.9cm (±1.27cm)
- Pelvic tilt (7°)
- Drain Pan: Stainless steel 1.5" deep (3.81cm)
- Large front storage drawer
- Side storage drawer – reversible

Patient Experience

- Removable, Seamless Upholstered Cushion
- Four lateral positions, adjustable length
- Comfort front leg cushion
- Heavy gauge steel extension leg pull-out to accommodate all patient heights

Electrical

- Hospital grade receptacle 110V with child safety cover (patient left)

Quality

- 3 year limited warranty
- Durable galvanized steel construction
- Easy-clean powder coat enamel surfaces
- International Models - 220/240V, 50/60Hz

COLORS

