



Hand
Control

Power Hi-Lo, Manual Back Exam Table with WheelBase[®], Programmable Hand Control, Foot Control
Model 4070-650-300

TABLE SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
Width: 28" (71.12cm)
Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
Width: 28" (71.12cm)
Length: 57" – 76" with foot section extended (144.7cm - 193cm)
Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
Stirrups: Four Lateral Positions, Adjustable Length
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60 HZ
International Model 4070-650-301: 220/240V, 50/60Hz
 4070-650-301: 220/240V, 50/60Hz

ACCESSORIES

042 Multi-Positioning Patient Assist Armrests / Safety Rails (Field Installable)
 251 Bierhoff Knee Crutches (Pair/Field installable)
 467 Mounting Plate for Welch Allen Light # 48955 (Patient Right)
 468 Mounting Plate for Welch Allen Light # 48955 (Patient Left)
 820 Stainless Steel Square Drain Pan (Factory installed only)
 822 Stainless Steel Drain Pan with Spout (Factory installed only)
 870 Drawer Warmer 110V (Factory Installed only)
 871 Hospital Grade Receptacle 110V (Factory installed only)
 900 Lifting bracket kit
 Fire retardant foam (CAL133) available

KEY FEATURES

- Power Hi-Lo with Air Spring Assisted Manual Backrest
- One-function Foot Control – Power Hi-Lo
- 650lb (294.84kg) Patient Weight Capacity
- OneTouch WheelBase[®] Integrated Central-Locking Caster System
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- Hand Control with OneTouch Patient Positioning[®] 2 Programmable Positions & Return-to-Home
- Removable Seamless Upholstered Cushion
- Side Storage Drawer
- Large Front Storage Drawer with Comfort Leg Pad
- Drawer Warmer
- Sturdy Steel Extension Leg Pullout with Seamless Ultra-Comfort Leg Cushion
- Stirrups
- Pelvic Tilt (7°)
- Stainless Steel Drain Pan
- Hospital Grade Electric Receptacle with Child Safety Cover (Patient Left)
- 3 Year Limited Warranty

COLORS

