

proglide 5016

podiatry/wound care procedure chair
Model: 5016-650-300



Shown with:
044 - Adjustable/Self Leveling Arms

COLORS



Procedure Chair, Power Leg Lift, Tilt, Hi-Lo & Back, Programmable Hand & Foot Controls

CHAIR SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat
Max. 37" at seat (93.9cm)
Seat Width: 28" (71.12cm)
Back Width: 26.5" (67.31cm)
Length: 73-80" with foot section extended
(185.42cm – 203.2cm)
Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60HZ - 220/240V, 50/60Hz
International Model (5016-650-301)

KEY FEATURES

- OneTouch WheelBase®
- OneTouch Patient Positioning® Programmable Foot Control for efficiencies
- OneTouch Patient Positioning® Programmable Hand Control
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Power Leg Lift makes it easier for patients, easier for staff
- Front Access to chair
- Ultra-Comfort Seat and Back Cushion
- Adjustable Articulating head rest provides ultimate comfort for patients
- 650 lb. Patient Weight Capacity
- Low access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- Hospital-Grade Electrical Receptacle 110v with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty
- Interchangeable Headrest Options



ACCESSORIES

044 - Adjustable/Self Leveling Arms (pair)
050 - Removable U-shaped Headrest
051 - Removable Rectangular Headrest
467 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Right
468 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Left

830 - Patient Safety Straps
900 - Lifting Bracket Kit
920 - Easy Installation Kit
Fire retardant foam (CAL133) available