

Manual Exam Table 5240





Manual Exam Table with Four Storage Drawers: Model 5240

SPECIFICATIONS

Height: 32.75" (83.2cm) Length: 58.5" - 74.5" (148.6cm - 189.23cm) Width: 28" (71.12cm) Stirrups: Infinitely adjustable Net Weight: 241lb (109kg) Shipping Weight: 320lb (145kg) Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

ACCESSORIES

- 010 Leg extension booster pad
- 037 Armboard assembly (specify patient right or left side, factory installed)
- 045 Accessory rail (pair, specify head or seat section, factory installed)
- 251 Bierhoff knee crutches (pair)
- 350 Pelvic tilt
- 365 I.V. pole (factory installed)
- 375 Safety rails assembly (pair, factory installed)
- 463 Light bracket capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- 701– Removable standard premium top
- 701F Removable standard premium top with CAL-133 foam
- 810 Restraint straps assembly (factory installed)
- 820 Stainless steel drain pan (factory installed)
- 822 Stainless steel drain pan with spout (factory installed)
- 849 Drawer warmer 220V (factory installed)
- 851 Drawer warmer 110V (factory installed)
- 852 Hospital grade receptacle 110V with child safe cover (factory installed, patient left side)
- 855 Hospital grade receptacle 220V with child safe cover (factory installed, patient left side)
- 5240F Fire retardant foam
- 5240B Base only

KEY FEATURES

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- Galvanized steel base, 18 gauge, powder coated
- Easy-reach articulating back adjustment
- One-touch, easy glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- Pull-out all steel non slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 2 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 5 year limited warranty

COLORS

