



Manual Exam Table with Four Storage Drawers: Model 5240

SPECIFICATIONS

Height: 32.75" (83.2cm)

Length: 58.5" – 74.5" (148.6cm – 189.23cm)

Width: 28" (71.12cm)

Stirrups: Infinitely adjustable

Net Weight: 241lb (109kg)

Shipping Weight: 320lb (145kg)

Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

ACCESSORIES

- 010 – Leg extension booster pad
- 037 – Armboard assembly (*specify patient right or left side, factory installed*)
- 045 – Accessory rail (*pair, specify head or seat section, factory installed*)
- 251 – Bierhoff knee crutches (*pair*)
- 350 – Pelvic tilt
- 365 – I.V. pole (*factory installed*)
- 375 – Safety rails assembly (*pair, factory installed*)
- 463 – Light bracket capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- 701 – Removable standard premium top
- 701F – Removable standard premium top with CAL-133 foam
- 810 – Restraint straps assembly (*factory installed*)
- 820 – Stainless steel drain pan (*factory installed*)
- 822 – Stainless steel drain pan with spout (*factory installed*)
- 849 – Drawer warmer 220V (*factory installed*)
- 851 – Drawer warmer 110V (*factory installed*)
- 852 – Hospital grade receptacle 110V with child safe cover (*factory installed, patient left side*)
- 855 – Hospital grade receptacle 220V with child safe cover (*factory installed, patient left side*)
- 5240F – Fire retardant foam
- 5240B – Base only

KEY FEATURES

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- Galvanized steel base, 18 gauge, powder coated
- Easy-reach articulating back adjustment
- One-touch, easy glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- Pull-out all steel non slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 2 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 5 year limited warranty

COLORS

