



Manual Exam Table with Five Drawers: Model 5250

## SPECIFICATIONS

- Height:** 32.75" (83.2cm)  
**Length:** 58.5" – 74.5" (148.6cm – 189.23cm)  
**Width:** 28" (71.12cm)  
**Stirrups:** Infinitely adjustable  
**Net Weight:** 241lb (109kg)  
**Shipping Weight:** 320lb (145kg)  
**Shipping Dimensions:** L64xW30xH41 (163cm x 76cm x 104cm)

## ACCESSORIES

- **010** – Leg extension booster pad
- **037** – Armboard assembly (*specify patient right or left side, factory installed*)
- **045** – Accessory rail (*pair, specify head or seat section, factory installed*)
- **251** – Bierhoff knee crutches (*pair*)
- **350** – Pelvic tilt
- **365** – I.V. pole (*factory installed*)
- **375** – Safety rails assembly (*pair, factory installed*)
- **463** – Light Bracket Capabilities (DOES NOT INCLUDE ACTUAL BRACKET)
- **701** – Removable standard premium top
- **701F** – Removable standard premium top with CAL-133 foam
- **810** – Restraint straps assembly (*factory installed*)
- **820** – Stainless steel drain pan (*factory installed*)
- **822** – Stainless steel drain pan with spout (*factory installed*)
- **849** – Drawer warmer 220V (*factory installed*)
- **851** – Drawer warmer 110V (*factory installed*)
- **852** – Hospital grade receptacle 110V with child safe cover (*factory installed, patient left side*)
- **855** – Hospital grade receptacle 220V with child safe cover (*factory installed, patient left side*)
- **5250F** – Fire retardant foam
- **5250B** – Base only

## KEY FEATURES

- 500lb capacity (226.8kg)
- Removable ultra-comfort seamless top
- Galvanized 18 gauge, powder coated steel base
- Easy-reach articulating back adjustment
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cm x 22.9cm)
- Pull out all steel non-slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 3 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cm x 39.4cm x 10.2cm)
- 5 year limited warranty

## COLORS

