



## Head Centric Procedure Chair with Programmable Hand and Foot Controls: Model 4010-650-200

### CHAIR SPECIFICATIONS

**sHeight:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat  
Max. 37" at seat (93.9cm)  
**Seat Width:** 28" (71.12cm)  
**Back Width:** 24.5" (62.23cm)  
**Length:** 70" – 80" with foot section extended (177.8cm – 203.2cm)  
**Leg Rest:** 28" W x 13" D (71.12cm W x 33.02cm D)  
**Dual Electro-Mechanical Actuator System** for Efficient Speed  
**Electrical:** 120 VAC, 60HZ  
**International Model 4010-650-201:** 220/240V, 50/60Hz

### KEY FEATURES

- Chair Designed for Head and Neck Centric Practices
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" ( ±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
- 2 Programmable Positions & Return-to-Home
- Foot Control with OneTouch Patient Positioning®
- 2 Programmable Positions and Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

### ACCESSORIES

- 044 Adjustable/self leveling arms pair- field installable
- 050 Removable U-shaped headrest
- 051 Removable rectangular headrest
- 467 Mounting plate for Welch Allyn Light #48955 - patient right
- 468 Mounting plate for Welch Allyn Light #48955 - patient left
- 830 Patient safety straps
- 900 Lifting bracket kit
- 920 Easy installation kit
- Fire retardant foam (CAL133)) available

### COLORS

