



Ultra Procedure Chair with Wheelbase®, Programmable Hand and Foot Controls: Model 4011-650-300

SPECIFICATIONS

Height: Min. 18.5" (±.5") 46.9cm (±1.27cm) at seat Max. 37" at seat (93.9cm)

Seat Width: 28" (71.12cm)

Back Width: 26.5" (67.31cm)

Length: 73"-80" with foot section extended (185.42cm – 203.2cm)

Leg Rest: 28" W x 11.5" D (71.12cm W x 29.21cm D)

Stirrups: Four Lateral Positions, Adjustable Length

Dual Electro-Mechanical Actuator System for Efficient Speed

Electrical: 120 VAC, 60HZ

International Model 4011-650-301: 220/240V, 50/60Hz

Net Weight: 495lb (224kg)

Shipping Weight: 535lb (242kg)

Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

ACCESSORIES

044 Adjustable / Self-Leveling Arms (Pair)

050 Removable U-Shaped Headrest

051 Removable Rectangular Headrest

052 Removable Oval Headrest

251 Bierhoff Knee Crutches (Pair)

467 Welch Allyn Light Bkt For #48955 (Patient Right)

468 Welch Allyn Light Bkt For #48955 (Patient Left)

830 Patient Safety Strap (Factory Installed Only)

873 220v Outlet Package (Factory Installed) (*International*)

874 220v Power Cord (All 650-101 Models) (*International*)

900 Lifting Bracket Kit

910 Easy Power Table Installation Kit

Fire Retardant Foam (Cal133) Available

KEY FEATURES

- OneTouch WheelBase® Integrated Central-Locking Caster System
- 650 lb. (294.84kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Foot Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Removable Seamless Upholstered Cushion
- Double-Articulating Headrest with Adjustable Height
- Interchangeable Headrest Options
- Stirrups
- Stainless Steel Drain Pan
- Trendelenburg, Lithotomy, Semi-Fowler Positions
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

COLORS

