



Stainless Steel One Step Foot Stool: Model SS8372

## CHAIR SPECIFICATIONS

- **Width: 30" (76.2cm)**
- **Height: 7.75" (19.7cm)**
- **Depth: 12" (30.5cm)**

## KEY FEATURES

- High quality, welded stainless steel construction
- 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- Non-slip rubber feet on stool legs
- Non-slip rubber step treads
- 3 year limited warranty

## ACCESSORIES

- N/A